

# TRANS-PHARMA™

## AUTOCLAVE BAGS & BIOHAZARD WASTE DISPOSAL

### AUTOCLAVE AND BIOHAZARD BAGS THAT CAN TAKE THE HEAT

TRANS-PHARMA™ autoclave cart liners can withstand the heat, making the transition from an autoclave to a shredder economical and virtually seamless. Our biohazard bags come in custom or standard sizes and the industry standard colors of orange, red, or yellow. As an option they can contain a heat strip to provide visual proof the contents have been sterilized and imprinted to increase safety, communication and proper disposal. No matter what the project calls for, TRANS-PHARMA™ products provide a single source solution.


#### CART LINERS

- Our cart liners have been tested to make the transition from autoclave to shredder virtually seamless
- Extend the life cycle of every cart
- Reduce wear and tear on your dumper by dumping the treated waste easily
- Enhance perception of clean/safe operation: Maintain carts free of debris
- Eliminate a potential source of nuisance odors
- Eliminate labor required to periodically remove waste debris from carts
- Consistent full capacity utilization—eliminating waste build-up inside carts
- Reduce occupational exposure, potential injuries related to cleaning waste from carts

#### BIO HAZARD BAGS—THE RIGHT SIZE FOR THE RIGHT JOB

Transcendia's TRANS-PHARMA™ bags come in a variety of sizes and pack options.

- Reduce cost by getting the right bag for the right purpose
- Get the right size! Bags that fit your containers—no more bags just a bit too small, or way too big
- Buy only the features you need on your bags
- Get your bags when you need them



#### APPLICATIONS:

Cart or Bin Liners  
Can Liners  
Medical Office Bags

Medical  
Food Services



# TRANS-PHARMA™

## AUTOCLAVE BAGS & BIOHAZARD WASTE DISPOSAL

### AUTOCLAVE BAGS & BIOHAZARD WASTE DISPOSAL

#### Procedure For Research And Teaching Laboratories

If your research or teaching activities generate waste that is autoclaved, please:

- Place all biohazardous waste to be autoclaved in a red or orange bag
- Prior to autoclaving, crisscross the bag's biohazard symbol and/or markings with heat sensitive autoclave tape
- Alternatively, mark the autoclaved bag with an "Autoclaved/Decontaminated" sticker prior to disposal in the normal trash. Because they are more visible, heat-sensitive tape or a sticker is required even if you use biohazardous waste bags that change to reveal the word "autoclaved" after autoclaving.
- Do not use red or orange bags for any other waste

### AVAILABILITY

Have your bags packaged the way you want.

- Individual Pouches



- By the Box



- On a roll

